

# RK Series Refrigerators

Rollex introduces the  
RK Series of full Stainless Steel Medical Refrigerators

The RK Series includes 2 sizes.

600 litre single door

739W x 805D x 2016H (mm)

1100 litre twin door

1260W x 805D x 2016H (mm)

Both designs come with either solid or glass doors. They have a beautifully finished stainless steel exterior and interior, giving them greater temperature retention.

Features:

- Digital display and controller
- Fan forced
- Self closing doors
- Stainless steel construction
- Easy to move with castors
- Lockable (solid door only)



RK600S: 600 Litre Solid Door



RK600G: 600 Litre Glass Door



