



The next generation blood managing system

SAFE AND RELIABLE REMOTE BLOOD MANAGEMENT

What is Hemosafe®

A new blood bank automation system to storage, preserve, distribute and deliver blood bags assuring full traceability and security.

- Easy to use (touch screen, vocal messages, kiosk features)
- Connected in real time with Blood Bank Information System
- Remote controllable from Blood Bank

Main advantages

- Automated distribution in remote locations
- Integrated access control and complete log of user activities
- Enhanced security provided by electronic crossmatch
- Reduces transportation costs
- Prevents blood expiration (FIFO)
- Global standards and international guidelines adopted (ISBT 128, AABB)

Technical data

- Load capacity: 150 bags (450 ml)
- Load per shelf: 15 bags (10 rotating shelves)
- Working temperature: +4°C with alarm and recording system
- Overall size: 1045x875x1805 mm



BLOOD
CENTER/HOSPITAL



OPERATING
ROOM



EMERGENCY
ROOM



REMOTE
HOSPITAL



Dati e caratteristiche tecniche possono essere variate senza preavviso
Appearance and specifications are subject to change without notice

ANGELANTONI Industrie SpA
BS biomedical division
Loc. Cimacole, 464
06056 Massa Martana (PG), Italy
Tel. (+39) 075 89551
Fax (+39) 075 8955312
e-mail: biomedical@angelantoni.it
www.angelantoni.it