

LABORATORY FREEZERS. The laboratory freezers from Lec Medical are available in four different sizes. From 50 Litres to 230 Litres with all four models having matching laboratory refrigerators available.

Control

The integral control panel on the outside of the fridge sits flush with the door and its easy to use design makes checking the internal temperature straightforward. As there is no need to open the door to check the temperature it keeps the internal air temperature more stable.

Security

The spark-free interior reduces the risk of internal explosion.

Security is essential in laboratories so all Lec Medical freezers come complete with door locks and two keys ensuring that only authorised individuals are able to access them.

Flexibility

Reversible doors across all models give greater flexibility where space is at a premium or where positioning the product is specific.

Adjustable feet are available on most models making it easier when positioning the product on uneven surfaces.

Environmentally Friendly

All our products are CFC free with many being HFC free too. Lec Medical refrigeration products are some of the most efficient products in the market using both efficient refrigerants and high specification insulation to help drive down the energy cost for our customers.

ISU27

SKU 444448572

Product dimensions (mm) 500w x 510d x 562h
 Packed dimensions (mm) 520w x 595d x 590h
 Product weight (kg) 17
 Packed weight (kg) 20
 Net capacity (litres) 38
 Gross capacity (litres) 50

ISU37

SKU 444449327

Product dimensions (mm) 495w x 580d x 830h
 Packed dimensions (mm) 560w x 655d x 850h
 Product weight (kg) 33
 Packed weight (kg) 36.5
 Net capacity (litres) 87
 Gross capacity (litres) 102

ISU57

SKU 444448522

Product dimensions (mm) 558w x 600d x 863h
 Packed dimensions (mm) 600w x 610d x 900h
 Product weight (kg) 38
 Packed weight (kg) 39.5
 Net capacity (litres) 94
 Gross capacity (litres) 116

ISU97

SKU 444448523

Product dimensions (mm) 595w x 600d x 1550h
 Packed dimensions (mm) 625w x 700d x 1610h
 Product weight (kg) 75
 Packed weight (kg) 77.5
 Net capacity (litres) 202
 Gross capacity (litres) 230

ISU27



ISU37



ISU57



ISU97

