

LABORATORY FRIDGE-FREEZERS. The laboratory fridge-freezers from Lec Medical are available in three different sizes, from 83.5 Litres to 246 Litres.

Control

The integral control panel on the outside of the fridge sits flush with the door and its easy to use design makes checking the internal temperature straightforward.

As there is no need to open the door to check the temperature it keeps the internal air temperature more stable.

Automatic defrost on the refrigerator area ensures more efficient cooling.

Security

The spark-free interior reduces the risk of internal explosion.

Security is essential in laboratories so all Lec Medical fridge-freezers come complete with door locks and two keys ensuring that only authorised individuals are able to access them.

Flexibility

Reversible doors across all models give greater flexibility where space is at a premium or where positioning the product is specific.

Adjustable feet are available on most models making it easier when positioning the product on uneven surfaces.

All of the laboratory fridge-freezers come with removable shelves and door furniture* giving maximum storage capacity.

Environmentally Friendly

All our products are CFC free with many being HFC free too. Lec Medical refrigeration products are some of the most efficient products in the market using both efficient refrigerants and high specification insulation to help drive down the energy cost for our customers.

ISR27

SKU	444448571
Product dimensions (mm)	502w x 540d x 660h
Packed dimensions (mm)	520w x 595d x 675h
Product weight (kg)	25
Packed weight (kg)	27
Fridge net capacity (litres)	75
Freezer net capacity (litres)	5.5
Fridge gross capacity (litres)	77
Freezer gross capacity (litres)	6.5

ISR57

SKU	444448518
Product dimensions (mm)	558w x 600d x 863h
Packed dimensions (mm)	600w x 610d x 900h
Product weight (kg)	36
Packed weight (kg)	37.5
Fridge net capacity (litres)	110
Freezer net capacity (litres)	17
Fridge gross capacity (litres)	112
Freezer gross capacity (litres)	18

IST47

SKU	444448520
Product dimensions (mm)	560w x 600d x 1572h
Packed dimensions (mm)	600w x 610d x 1620h
Product weight (kg)	68
Packed weight (kg)	70
Fridge net capacity (litres)	152
Freezer net capacity (litres)	79
Fridge gross capacity (litres)	157
Freezer gross capacity (litres)	89



ISR27



ISR57



IST47