

ROLLEX MEDICAL (AUST) PTY LTD
 PO Box 170 | Mt Kuring-gai 2080 | Sydney, NSW | Australia
 Phone: +61 2 9457 6144 | Fax: +61 2 9457 6516
 Email: sales@rollexmedical.com | www.rollexmedical.com.au

Model Variations

MATOS Aria Eco 374 R – Solid Door
 MATOS Aria Eco 374 R/G – Glass Door



MATOS Aria Eco 374 R/G

The MATOS® Aria Eco fridge series is a new range of cost-effective European-built Medical and Vaccination Refrigerators featuring the latest technology in controllers and compressors for accurate temperature stability.

Advanced Controller

Featuring the Carel soft touch controller - accurate to 0.1°C with min/max readings at the touch of a button.

Lock Fitted as Standard

Maintains security and ensures the unit can only be opened by authorised personnel.

Cost-Effective

Built with a High Impact Polystyrene Interior to provide a robust and in-expensive range of refrigerators.

Solid or Glass Door Options

Solid or glass door options enabling you to select the format which best suits your requirements. Glass doors are double glazed, toughened safety-glass with argon gas between panels for better insulation.

Auto Defrost

The advanced controller initialises a regular defrost sequence to decrease the chance of ice build-up while maintaining the temperature below 8°C.

Extended Warranty

Standard 12 month warranty however an extended warranty option is available subject to Terms and Conditions

CCTC Certified

CCTC and QCCP Approved. Exclusively endorsed by The Pharmacy Guild of Australia.

Endorsed by



The Pharmacy
 Guild of Australia

Specifications

Capacity	342L Net / 358L Gross	Power usage	Solid 4.3 kW Glass 4.8 kW / 24 hours
Dimensions external	600w x 640d x 2040h	Air convection	Fan assisted
Dimensions internal	505w x 462d x 1555h	Temperature range	1°C to +10°C
Door	Solid and Glass available Reversible	Controller	Microprocessor with external LCD graphic display
Weight	Solid 75kg / Glass 85kg	Interior/Exterior	(HIPS) High impact polystyrene
Shelves fitted	6	Temperature resolution	0.1°C
Castors	2 wheels rear, 2 adjustable feet front	Temperature fluctuation at 4°C	+/- 0.5°C
Internal light	Yes	Refrigerant	R134A
Door Open Alarm	Yes	Power	230V / 50Hz
Noise Level dB	46.8		