

ROLLEX MEDICAL (AUST) PTY LTD  
 PO Box 170 | Mt Kuring-gai 2080 | Sydney, NSW | Australia  
 Phone: 1300 880 441 | Fax: 1330 880 451  
 Email: sales@rollexmedical.com | www.rollexmedical.com.au

The Matos® Plus Eco freezer series is a new range of purpose built Medical Freezers utilising patented technology. The cabinet and componentry have been specially designed for use in the Medical field.

## Model Variations

Eco 85 F – Solid Door



MATOS Eco 85 F

## Interior Finish

All Matos® Plus Eco freezers come with a high quality stainless steel interior and use polyurethane foam thermal insulation.

## Solid or Glass Door Options

Most models come with solid or glass door options enabling you to select the format which best suits your requirements.

## Monitoring capability

Continuous temperature monitoring with sms and email alerts available – ensures a robust complete monitoring record of your product – even with loss of power.

## Lock Fitted as Standard

Maintains security and ensures the unit can only be opened by authorised personnel.

## Energy Efficient

Matos® Plus Eco freezers are more energy efficient and kinder to the environment, meaning the running costs of the freezer are kept as low as possible.

## 3 Year Extended Warranty

When serviced annually by an approved service technician,

## Specifications

Capacity	85L	Air convection	Natural
Dimensions External	620w x 880h x 650d	Temperature range	-25°C to 0°C
Dimensions Internal	420w x 590h x 395d	Controller	Microprocessor with external LCD graphic display
Max shelf workload	10kg	Interior	Stainless Steel
Max unit workload	30kg	Housing	Powder coated sheet (s/steel optional)
Door	Solid	Colour	Silver/Grey
Weight	60kg	Temperature accuracy	1°C
Shelves fit/max	2/4	Temperature fluctuation at -10°C	+/- 0.4°C
Refrigerant	R290	Voltage 50/60Hz	115/230V
Refrigerant amount	0.080	Internal light	No
Castors	No	Door Open Alarm	Yes
Power	200W	Access Port	Yes