

ROLLEX MEDICAL (AUST) PTY LTD
 PO Box 170 | Mt Kuring-gai 2080 | Sydney, NSW | Australia
 Phone: +61 2 9457 6144 | Fax: +61 2 9457 6516
 Email: sales@rollexmedical.com | www.rollexmedical.com.au

Model Variations

MATOS Aria Eco 374 R – Solid Door
 MATOS Aria Eco 374 R/G – Glass Door



MATOS Aria Eco 374 R/G

The MATOS® Aria Eco fridge series is a new range of cost-effective European-built Medical and Vaccination Refrigerators featuring the latest technology in controllers and compressors for accurate temperature stability.

Advanced Controller

Featuring the Carel soft touch controller - accurate to 0.1°C with min/max readings at the touch of a button.

Lock Fitted as Standard

Maintains security and ensures the unit can only be opened by authorised personnel.

Cost-Effective

Built with a High Impact Polystyrene Interior to provide a robust and in-expensive range of refrigerators.

Solid or Glass Door Options

Solid or glass door options enabling you to select the format which best suits your requirements. Glass doors are double glazed, toughened safety-glass with argon gas between panels for better insulation.

Auto Defrost

The advanced controller initialises a regular defrost sequence to decrease the chance of ice build-up while maintaining the temperature below 8°C.

Extended Warranty

Standard 12 month warranty however an extended warranty option is available subject to Terms and Conditions

CCTC / QCPP Certified

CCTC and QCPP Approved Vaccine Pharmacy Graded Refrigerator.



Specifications

| | | | |
|---------------------|---|--------------------------------|---|
| Capacity | 345L Net / 372L Gross | Power usage | Solid 4.3 kW Glass 4.8 kW / 24 hours |
| Dimensions external | 597w x 645d x 1896h | Air convection | Fan assisted |
| Dimensions internal | 505w x 462d x 1555h | Temperature range | 1°C to +10°C |
| Door | Solid and Glass available Reversible | Controller | Microprocessor with external LCD graphic display |
| Weight | Solid 70kg / Glass 80kg | Interior/Exterior | (HIPS) High impact polystyrene |
| Shelves fitted | 6 | Temperature resolution | 0.1°C |
| Castors | 2 wheels rear, 2 adjustable feet front | Temperature fluctuation at 4°C | +/- 0.5°C |
| Internal light | Yes | Refrigerant | R134A |
| Door Open Alarm | Yes | Power | 230V / 50Hz |
| Noise Level dB | 46.8 | | |