

ROLLEX MEDICAL (AUST) PTY LTD
 PO Box 170 | Mt Kuring-gai 2080 | Sydney, NSW | Australia
 Phone: 1300 880 441 | Fax: 1330 880 451
 Email: sales@rollexmedical.com | www.rollexmedical.com.au

The Matos® Eco freezer series is a new range of purpose built Medical Freezers utilising patented technology. The cabinet and componentry have been specially designed for use in the Medical field.

Model Variations

Eco 125 F – Solid Door



MATOS Eco 125 F

Interior Finish

All Matos® Eco freezers come with a high quality stainless steel interior and use polyurethane foam thermal insulation.

Solid or Glass Door Options

Most models come with solid or glass door options enabling you to select the format which best suits your requirements.

Monitoring capability

Continuous temperature monitoring with sms and email alerts available – ensures a robust complete monitoring record of your product – even with loss of power:

Lock Fitted as Standard

Maintains security and ensures the unit can only be opened by authorised personnel.

Energy Efficient

Matos® Eco freezers are more energy efficient and kinder to the environment, meaning the running costs of the freezer are kept as low as possible.

3 Year Extended Warranty

When serviced annually by an approved service technician, the warranty is extended to 3 years.

Specifications

Capacity	125L	Air convection	Natural
Dimensions External	675w × 1170h × 800d	Temperature range	-40°C to 0°C
Dimensions Internal	396w × 600h × 524d	Controller	Microprocessor with external LCD graphic display
Max shelf workload	10kg	Interior	Stainless Steel
Max unit workload	50kg	Housing	Powder coated sheet (s/steel optional)
Door	Solid	Colour	Silver/Grey
Weight	90kg	Temperature accuracy	1°C
Shelves fit/max	2/3	Temperature fluctuation at -10°C	+/- 0.7°C
Refrigerant	R290	Voltage 50/60Hz	115/230V
Refrigerant amount	0.18kg	Internal light	No
Castors	No	Door Open Alarm	Yes
Power	450W	Access Port	Yes