

ROLLEX MEDICAL (AUST) PTY LTD  
 PO Box 170 | Mt Kuring-gai 2080 | Sydney, NSW | Australia  
 Phone: 1300 880 441 | Fax: 1330 880 451  
 Email: sales@rollexmedical.com | www.rollexmedical.com.au

## Model Variations

- Eco 1460 R/DT – Solid Door, hot gas defrost
- Eco 1460 R/GDT – Glass Door, hot gas defrost



MATOS PLUS Eco 1460 R/GDT



## Specifications

|                     |                      |                                |  |
|---------------------|----------------------|--------------------------------|--|
| Capacity            | 1460L                | Power                          | 550W   |
| Dimensions External | 1450w x 1970h x 950d | Air convection                 | Forced Air                                       |
| Dimensions Internal | 1340w x 1460h x 750d | Temperature range              | 0°C – +10°C                                      |
| Max shelf workload  | 30kg                 | Controller                     | Microprocessor with external LCD graphic display |
| Max unit workload   | 300kg                | Interior                       | Aluminium (s/steel optional)                     |
| Door                | Solid/Glass option   | Housing                        | Powder coated sheet (s/steel optional)           |
| Weight              | 200kg                | Temperature resolution         | 0.1°C  |
| Shelves fit/max     | 2 x 3/11             | Temperature fluctuation at 4°C | +/- 0.5°C  |
| Gas defrost         | Optional             | Voltage 50/60Hz                | 115/230V   |
| Refrigerant         | R404A                | Internal light                 | Yes  |
| Refrigerant amount  | 0.27kg               | Door Open Alarm                | Yes  |
| Castors             | Yes                  | Access port                    | Yes  |

The Matos® Plus Eco fridge series is a new range of purpose built Medical and Vaccination Refrigerators utilising patented technology. The cabinet and componentry have been specially designed for use in the Medical field.

### Interior Finish

All Matos® Plus Eco fridges come with a high quality aluminium interior and use polyurethane foam thermal insulation. A stainless steel interior is also offered as an option.

### Lock Fitted as Standard

Maintains security and ensures the unit can only be opened by authorised personnel.

### Solid or Glass Door Options

Most models come with solid or glass door options enabling you to select the format which best suits your requirements.

### Auto Defrost

The advanced controller initializes a regular defrost sequence to decrease the chance of ice build-up while maintaining the temperature below 8°C. Matos® Plus Eco fridges can also be ordered with a gas defrost capability.

### Energy Efficient

Matos® Plus Eco fridges are more energy efficient and kinder to the environment, meaning the running costs of the refrigerator are kept as low as possible.

### 3 Year Extended Warranty

When serviced by an approved service technician, the warranty is extended to 3 years.

### CCTC / QCPP Certified

CCTC and QCPP Approved Vaccine Pharmacy Graded Refrigerator.