

ROLLEX MEDICAL (AUST) PTY LTD  
 PO Box 170 | Mt Kuring-gai 2080 | Sydney, NSW | Australia  
 Phone: 1300 880 441 | Fax: 1330 880 451  
 Email: sales@rollexmedical.com | www.rollexmedical.com.au

*The Matos® Plus Eco Cooled Incubators are a new range of purpose built Medical Chambers utilising patented technology. The cabinet and componentry have been specially designed for use in the Medical field.*

## Model Variations

- Eco 68 S – Solid Door
- Eco 68 S/G – Glass Door



MATOS Eco 68 S



MATOS Eco 68 S/G

## Interior Finish

All Matos® Plus Eco **Cooled Incubators** come with a high quality aluminium interior and use polyurethane foam thermal insulation. A stainless steel interior is also offered as an option.

## Lock Fitted as Standard

Maintains security and ensures the unit can only be opened by authorised personnel.

## Solid or Glass Door Options

Most models come with solid or glass door options enabling you to select the format which best suits your requirements.

## Auto Defrost

The advanced controller initializes a regular defrost sequence to decrease the chance of ice build-up while maintaining the temperature.

## Energy Efficient

Matos® Plus Eco **Cooled Incubators** are more energy efficient and kinder to the environment, meaning the running costs of the refrigerator are kept as low as possible.

## 3 Year Extended Warranty

When serviced by an approved service technician, the warranty is extended to 3 years.

## Specifications

Capacity	68L	Power	160W
Dimensions External	570w x 600h x 680d	Air convection	Forced Air
Dimensions Internal	470w x 430h x 300d	Temperature range	+3°C – +40°C
Max shelf workload	10kg	Controller	Microprocessor with external LCD graphic display
Max unit workload	20kg	Interior	Aluminium (s/steel optional)
Door	Solid/Glass option	Housing	Powder coated sheet (s/steel optional)
Weight	30kg	Temperature resolution	0.1°C
Shelves fit/max	2/2	Temperature fluctuation at 4°C	+/- 0.3°C
Gas defrost	Optional	Voltage 50/60Hz	115/230V
Refrigerant	R1234	Door Open Alarm	Yes
Refrigerant amount	55 grams	Access Port	Yes
Castors	Optional	Internal light	Yes