




COUNTERTOP LABORATORY FREEZER SOLID DOOR 50L WHITE



- ✓ **-18°C to -20°C Operating Temperature Range**
Conforms to MHRA guidelines
- ✓ **NEW Unique Anti Microbial Handle**
Helps prevent the growth of bacteria and minimise the potential for cross contamination
- ✓ **Fitted Lock as Standard**
Maintains security and ensures the unit can only be accessed by approved personnel
- ✓ **Wire Shelf**
Enables optimum air flow inside unit



The Pharmacy
Guild of Australia

SKU CODE:

PRODUCT CODE:

DESCRIPTION:

444442352

ISU27CAU

Countertop Laboratory Freezer Solid Door 50L White



FULL SPECIFICATION – ISU27CAU COUNTERTOP LABORATORY FREEZER SOLID DOOR 50L WHITE

Description	Countertop
SKU Code	444442352
EAN Code	5052263023529
GDPA Code	ISU27CAU
Gross Capacity (L)	50
Net Capacity (L)	38
Unit External Width (mm)	500
Unit External Depth (mm)	510
Unit External Height (mm)	562
Product Weight (kg)	17
Packed Width (mm)	520
Packed Depth (mm)	595
Packed Height (mm)	590
Packed Weight (kg)	20
Country of Origin	UK
Internal Dimensions (mm) (WxDxH)	345 x 330 x 365
External Finish	White
Internal Finish	White
Type of Plastic	SRL800 (HIPS/ABS Mix)
Temperature Range	-18°C to -20°C
Colour	White/Green
Door Type	Solid
Door Style	Hinged, Reversible
Door Open Depth (mm)	976
Replacement Door Seal	No
Open Door Sensor	No
Lock	Yes
Adjustable Feet	No
Temperature Display on Panel	Yes
Temperature Alarm	Yes
Number of Shelves	1
Number of Shelf Positions	Fixed
Shelf Size	-
Shelf Material	Wire
Moulded Drawers	0
Adjustable Door Bins	0
Interior Light	No
Fan	No
Spark Free	Yes
Auto or Manual defrost	Manual
Controller type	Digital, External
Climate Class	N
Refrigerant Type	R600a
CFC Free	Yes
HFC Free	No
HCFC Free	Yes
Australian Standard Plug	Yes
Other Plug Types Available	Yes
Voltage Supply	230
Aux Output	Yes
Insulation Thickness (mm)	30
Insulation Type	Cyclopentane
Anti Microbial Handle	Yes
Tariff Code	8418402090