



UNDER-COUNTER LABORATORY FREEZER SOLID DOOR 151L WHITE



- ✓ **-15°C to -18°C Operating Temperature Range**
- ✓ **NEW Advanced Controller**
Records temperatures accurate to 0.1°C
- ✓ **NEW Unique Anti Microbial Handle**
Helps prevent the growth of bacteria and minimise the potential for cross contamination
- ✓ **Fitted Lock as Standard**
Maintains security and ensures the unit can only be accessed by approved personnel
- ✓ **Fitted Rollers**
For ease of movement
- ✓ **Wire Baskets**
Enables optimum air flow inside the unit
- ✓ **Audible Warning Alarms**
Alerts user should internal temperature deviate from set parameters



The Pharmacy
Guild of Australia

| SKU CODE: | PRODUCT CODE: | DESCRIPTION: |
|-----------|---------------|--|
| 444440968 | LSF151 | Under-counter Laboratory Freezer Solid Door 151L White |



FULL SPECIFICATION – LSF151 UNDER-COUNTER LABORATORY FREEZER SOLID DOOR 151L WHITE



| | |
|---|-------------------------------------|
| Description | Under-counter |
| SKU Code | 444440968 |
| EAN Code | 5052263009684 |
| GDPA Code | LSF151 |
| Gross Capacity (L) | 151 |
| Net Capacity (L) | 134 |
| Unit External Width (mm) | 595 |
| Unit External Depth (mm) | 595 |
| Unit External Height (mm) | 845 |
| Product Weight (kg) | 36 |
| Packed Width (mm) | 657 |
| Packed Depth (mm) | 694 |
| Packed Height (mm) | 884 |
| Packed Weight (kg) | 40 |
| Country of Origin | UK |
| Internal Dimensions (mm) (WxDxH) | 500 x 450 x 705 (step: 165d x 185h) |
| External Finish | White |
| Internal Finish | White |
| Type of Plastic | SRL800 (HIPS/ABS Mix) |
| Temperature Range | -15°C to -18°C |
| Colour | White/Green |
| Door Type | Solid |
| Door Style | Hinged, Reversible |
| Door Open Depth (mm) | 1165 |
| Replacement Door Seal | No |
| Open Door Sensor | Yes |
| Lock | Yes |
| Adjustable Feet | 2 Adjustable Feet, 2 Rear Rollers |
| Temperature Display on Panel | Yes |
| Temperature Alarm | Yes |
| Number of Shelves | 3 Baskets |
| Number of Shelf Positions | Fixed |
| Shelf Size | - |
| Shelf Material | Wire Baskets |
| Moulded Drawers | 0 |
| Adjustable Door Bins | 0 |
| Interior Light | No |
| Fan | No |
| Spark Free | Yes |
| Auto or Manual defrost | Manual |
| Controller type | Digital, External |
| Climate Class | SN |
| Refrigerant Type | R600a |
| CFC Free | Yes |
| HFC Free | Yes |
| HCFC Free | Yes |
| UK 3 Pin Plug | Yes |
| Other Plug Types Available | Yes |
| Voltage Supply | 230 |
| Aux Output | Yes |
| Insulation Thickness (mm) | 45 |
| Insulation Type | Cyclopentane |
| Anti Microbial Handle | Yes |