



CLT383 LOW TEMPERATURE CHEST FREEZER

-10°C to -45°C



- ✓ **-10°C to -45°C set temperature**
Ideal for medium term sample storage.
- ✓ **External digital controller**
For temperature readings without having to open the door, so preventing unnecessary temperature fluctuations within the freezer.
- ✓ **Lock and two keys**
For improved safety, accountability and security over domestic refrigeration, a lockable fridge helps prevent unauthorised access to valuable and dangerous contents.
- ✓ **Spark free**
There is nothing electrical inside the freezer so it is safe to store reactive substances without the risk of explosion.
- ✓ **Audible and visual high and low temperature alarms**
To alert anybody nearby as soon as the temperature begins to go too high or low, so that action can be taken to rectify the problem and prevent damage to the contents.



www.facebook.com/Lec-Medical
www.twitter.com/Lec_Medical

SKU CODE:	PRODUCT CODE:	DESCRIPTION:
444440162	CLT383	Low Temperature Chest Freezer -10°C to -45°C



FULL SPECIFICATION – CLT383 Low Temperature Chest Freezer -10°C to -45°C



SKU Code	444440162
GDPA Code	CLT383
Gross Capacity (L)	383
Unit External Width (mm)	1560
Unit External Depth (mm)	600
Unit External Height (mm)	855
Packed Width (mm)	1600
Packed Depth (mm)	700
Packed Height (mm)	960
Internal Dimensions (mm)	1400 x 440 x 650
Product Weight (kg)	80
Packed Weight (kg)	95
Colour	White
Temperature Range	-10°C to -45°C
WEEE Class	Medical Devices - Freezers
Style	Chest
Door Style	Solid, Hinged Lid
Replaceable Door Seal	No
Open Door Sensor	No
Lock	Yes
Adjustable Feet	Fitted Castors
Temperature Display	Yes
Audible and Visual Temperature Alarm	Yes
Noise Level	N/A
Shelves/Baskets	1 basket
Interior Light	No
Power Consumption (kw/24h)	5.2
Control	Electronic
Climate Class	N
Refrigerant Type	R414
AUX output	Yes
UK 3-pin plug	Yes
Voltage Supply	220-240v / 50Hz
Spark Free	Yes