



## FREESTANDING LABORATORY FREEZER SOLID DOOR 232L WHITE



- ✓ **-18°C to -25°C Operating Temperature Range**  
Desired temperature can be set by the user
- ✓ **NEW Advanced Controller**  
Records temperatures accurate to 0.1°C
- ✓ **NEW Unique Anti Microbial Handle**  
Helps prevent the growth of bacteria and minimise the potential for cross contamination
- ✓ **Fitted Lock as Standard**  
Maintains security and ensures the unit can only be accessed by approved personnel
- ✓ **Fitted Rollers**  
For ease of movement
- ✓ **Plastic Drawers and Flaps**  
Maximises space and offers protection to contents
- ✓ **Audible Warning Alarms**  
Alerts user should internal temperature deviate from set parameters



[www.facebook.com/Lec-Medical](https://www.facebook.com/Lec-Medical)  
[www.twitter.com/Lec\\_Medical](https://www.twitter.com/Lec_Medical)

SKU CODE:

PRODUCT CODE:

DESCRIPTION:

444440970

LSF232UK

Freestanding Laboratory Freezer Solid Door 232L White



## FULL SPECIFICATION – LSF232UK FREESTANDING LABORATORY FREEZER SOLID DOOR 232L WHITE



<b>Description</b>	Freestanding
<b>SKU Code</b>	444440970
<b>EAN Code</b>	5052263009707
<b>GDPA Code</b>	LSF232UK
<b>Gross Capacity (L)</b>	232
<b>Net Capacity (L)</b>	212
<b>Unit External Width (mm)</b>	595
<b>Unit External Depth (mm)</b>	634
<b>Unit External Height (mm)</b>	1570
<b>Product Weight (kg)</b>	60
<b>Packed Width (mm)</b>	657
<b>Packed Depth (mm)</b>	734
<b>Packed Height (mm)</b>	1585
<b>Packed Weight (kg)</b>	65
<b>Country of Origin</b>	UK
<b>Internal Dimensions (mm) (WxDxH)</b>	432 x 430 x 1300 (step: 190d x 255h)
<b>External Finish</b>	White
<b>Internal Finish</b>	White
<b>Type of Plastic</b>	SRL800 (HIPS/ABS Mix)
<b>Temperature Range</b>	-18°C to -25°C
<b>Colour</b>	White/Green
<b>Door Type</b>	Solid
<b>Door Style</b>	Hinged, Reversible
<b>Door Open Depth (mm)</b>	1185
<b>Replacement Door Seal</b>	No
<b>Open Door Sensor</b>	Yes
<b>Lock</b>	Yes
<b>Adjustable Feet</b>	2 Adjustable Feet, 2 Rear Rollers
<b>Temperature Display on Panel</b>	Yes
<b>Temperature Alarm</b>	Yes
<b>Number of Shelves</b>	4 Drawers, 2 Fixed Flaps
<b>Number of Shelf Positions</b>	Fixed
<b>Shelf Size</b>	-
<b>Shelf Material</b>	Wire Baskets
<b>Moulded Drawers</b>	4 Solid
<b>Adjustable Door Bins</b>	0
<b>Interior Light</b>	No
<b>Fan</b>	No
<b>Spark Free</b>	Yes
<b>Auto or Manual defrost</b>	Manual
<b>Controller type</b>	Digital, External
<b>Climate Class</b>	SN
<b>Refrigerant Type</b>	R600a
<b>CFC Free</b>	Yes
<b>HFC Free</b>	Yes
<b>HCFC Free</b>	Yes
<b>UK 3 Pin Plug</b>	Yes
<b>Other Plug Types Available</b>	Yes
<b>Voltage Supply</b>	230
<b>Aux Output</b>	Yes
<b>Insulation Thickness (mm)</b>	60
<b>Insulation Type</b>	Cyclopentane
<b>Anti Microbial Handle</b>	Yes