







FREESTANDING PHARMACY FRIDGE SOLID DOOR 444L WHITE



- 
+2°C to +8°C Operating Temperature Range
 Conforms to MHRA guidelines
- 
Fitted Lock as Standard
 Maintains security and ensures the unit can only be accessed by approved personnel
- 
Fitted Castors
 For ease of movement
- 
Wire Shelves
 Enables optimum air flow inside the unit



The Pharmacy
Guild of Australia

SKU CODE:

PRODUCT CODE:

DESCRIPTION:

444442342

PE1607CAU

Freestanding Pharmacy Fridge Solid Door 444L White



FULL SPECIFICATION – PE1607CAU FREESTANDING PHARMACY FRIDGE SOLID DOOR 444L WHITE

Description	Freestanding
SKU Code	444442342
EAN Code	5052263023420
GDPA Code	PE1607CAU
Gross Capacity (litres)	444
Net Capacity (litres)	393
Unit External Width (mm)	685
Unit External Depth (mm)	910
Unit External Height (mm)	1688
Product Weight (kg)	96
Packed Width (mm)	770
Packed Depth (mm)	950
Packed Height (mm)	1810
Packed Weight (kg)	117
Country of Origin	UK
Internal Dimensions (mm) (WxDxH)	590 x 670 x 1170
Temperature Range	+2°C to +8°C
Colour	White/Green
External Finish	Painted Steel
Internal Finish	HIPS*
Type of plastic	SRL800 (HIPS/ABS Mix)
Door Finish	Painted Steel
Door Type	Solid
Door Style	Hinged, R/H only
Door Open Depth (mm)	TBC
Replacement Door Seal	Yes
Open Door Sensor	Yes
Lock	Yes
Adjustable Feet (number & position)	Castors
Temperature Display on Panel	Yes
Temperature Alarm	Yes
Number of Shelves	5
Number of Shelf Positions	30
Shelf Size (mm) (WxDxH)	5: 536 x 609 x 6
Shelf Material	Wire
Number of Baskets	0
Basket Material	-
Interior Light	No
Fan	Yes
Spark Free	No
Auto or Manual Defrost	Automatic
Controller type	Digital, External
Climate Class	N
Refrigerant Type	R134a
CFC Free	Yes
HFC Free	No
HCFC Free	Yes
Ammonia Free	Yes
Australian Standard Plug	Yes
Other Plug Types Available	Yes
Voltage Supply	230
Aux Output	Yes
Insulation Thickness (mm)	40
Insulation Type	Cyclopentane
Anti Microbial Handle	No
Tariff Code	8418211000