



## FREESTANDING PHARMACY FRIDGE SOLID DOOR 273L WHITE



- 
**+2°C to +8°C Operating Temperature Range**  
 Conforms to MHRA guidelines
- 
**NEW Advanced Controller**  
 Records temperatures accurate to 0.1°C
- 
**NEW Unique Anti Microbial Handle**  
 Helps prevent the growth of bacteria and minimise the potential for cross contamination
- 
**Fitted Lock as Standard**  
 Maintains security and ensures the unit can only be accessed by approved personnel
- 
**Internal Light**
- 
**Fitted Rollers**  
 For ease of movement
- 
**Wire Shelves and Basket**  
 Enables optimum air flow inside the unit
- 
**Audible Warning Alarms**  
 Alerts user should internal temperature deviate from set parameters



The Pharmacy  
Guild of Australia

SKU CODE:	PRODUCT CODE:	DESCRIPTION:
444440952	PSR273AU	Freestanding Pharmacy Fridge Solid Door 273L White



**FULL SPECIFICATION – PSR273AU FREESTANDING PHARMACY FRIDGE SOLID DOOR  
273L WHITE**



<b>Description</b>	Freestanding
<b>SKU Code</b>	444440952
<b>EAN Code</b>	5052263009509
<b>GDPA Code</b>	PSR273AU
<b>Gross Capacity (litres)</b>	273
<b>Net Capacity (litres)</b>	269
<b>Unit External Width (mm)</b>	595
<b>Unit External Depth (mm)</b>	640 (w/o handle), 670 (with handle)
<b>Unit External Height (mm)</b>	1565
<b>Product Weight (kg)</b>	60
<b>Packed Width (mm)</b>	668
<b>Packed Depth (mm)</b>	737
<b>Packed Height (mm)</b>	1633
<b>Packed Weight (kg)</b>	65
<b>Country of Origin</b>	UK
<b>Internal Dimensions (mm) (WxDxH)</b>	480 x 445 x 1382 (step: 175d x 255h)
<b>Temperature Range</b>	+2°C to +8°C
<b>Colour</b>	White/Green
<b>External Finish</b>	Painted Steel
<b>Internal Finish</b>	HIPS*
<b>Type of plastic</b>	SRL800 (HIPS/ABS Mix)
<b>Door Finish</b>	Painted Steel
<b>Door Type</b>	Solid
<b>Door Style</b>	Right Hand Hinge as Default, Is reversible at an additional cost
<b>Door Open Depth (mm)</b>	1185
<b>Replacement Door Seal</b>	No
<b>Open Door Sensor</b>	Yes
<b>Lock</b>	Yes, 2 Keys
<b>Adjustable Feet (number &amp; position)</b>	2 Adjustable Front Feet, 2 Rear Rollers Yes
<b>Temperature Display on Panel</b>	Yes
<b>Temperature Alarm</b>	5
<b>Number of Shelves</b>	12
<b>Number of Shelf Positions</b>	4: 385 x 500 x 6,
<b>Shelf Size (mm) (WxDxH)</b>	1: 340 x 497 x 6
<b>Shelf Material</b>	Wire
<b>Number of Baskets</b>	1
<b>Basket Material</b>	Wire
<b>Interior Light</b>	Yes
<b>Fan</b>	Yes
<b>Spark Free</b>	No
<b>Auto or Manual Defrost</b>	Automatic
<b>Controller type</b>	Digital, External
<b>Climate Class</b>	SN
<b>Refrigerant Type</b>	R600a
<b>CFC Free</b>	Yes
<b>HFC Free</b>	Yes
<b>HCFC Free</b>	Yes
<b>Ammonia Free</b>	Yes
<b>Australian Standard Plug</b>	Yes
<b>Other Plug Types Available</b>	Yes
<b>Voltage Supply</b>	220-240v/50hz
<b>Aux Output</b>	Yes
<b>Insulation Thickness (mm)</b>	60
<b>Insulation Type</b>	Cyclopentane
<b>Anti Microbial Handle</b>	Yes
<b>Optional Drawer System Available</b>	Yes